

## Fruit Fly Exercise 7 – Level 1

### Goal

In this exercise, you will use StarGenetics, a software tool that simulates mating experiments, to analyze the nature and mode of inheritance of specific genetic traits.

### Learning Objectives

After completing this exercise, you will be able to:

1. Identify the phenotype of a fruit fly within StarGenetics.
2. Determine whether mutations are in the same or different genes through the analysis of results from genetic crosses.
3. Determine whether a phenotype is dominant or recessive relative to another phenotype through the analysis of results from genetic crosses.
4. Infer and assign genotypes of individual organisms using proper nomenclature of alleles.
5. Generate a genetic map that specifies the map distances between alleles.

### Getting started with StarGenetics

- To get to StarGenetics, please navigate to: <http://star.mit.edu/genetics/>.
- Click on the **Start** button to launch the application.
- Click **Trust** when a prompt appears asking if you trust the certificate.
- Click on **File → New** on the main menu.
- Click on the **Fruit Fly Exercise 7 – Level 1** file.

Fruit flies rely on their olfactory system to locate food, to find a mate, and to find the right place to lay their eggs. In fact, although fruit flies have brains smaller than a poppy seed, they are able to incorporate various types of visual and olfactory information to find food sources in a vast array of environments.

Fruit flies are able to differentiate among many different types of odors (they can be trained to avoid specific odors too!). In *Drosophila*, each antenna contains approximately 1200 olfactory neurons that help with the task of detecting and differentiating between various scents.

Fascinated by this fact, you decide to study the genetics of odor differentiation in fruit flies. You start with a fruit fly strain that has unusually large antennae. These flies were found on tropical islands in the Pacific and are known to feed on a wide range of fruits. These flies also avoid certain fruits that contain substances that are toxic to the flies. You hypothesize that the large antennae of these fruit flies might confer to them the ability to discriminate between the numerous edible and toxic fruits found on the Pacific islands.

To identify the genes that are important to large antennae development, you mutagenize wild-type flies to generate mutants with antennae defects. Via this process, you identify flies that have either shortened antennae or no antennae at all. In your **Strains** box in StarGenetics, you will find male (M) and female (F) flies for four mutant strains: **Mutant 1**, **Mutant 2**, **Mutant 3** and **Mutant 4**. In preliminary experiments, you observe that these mutant flies have a compromised “sense of smell”! You also find that each mutant fly contains a single mutation and is true breeding. You now need to perform further genetic analysis of these mutants.

- 1 Describe the phenotypes of Mutant 1, Mutant 2, Mutant 3 and Mutant 4.

**Answer**

**2** Are the mutations in Mutant 1, Mutant 2, Mutant 3, and Mutant 4 in the same gene or in different genes? Describe the reason for all crosses performed to answer this question, their outcomes, and the conclusion that you drew from each of these crosses.

- You can set up a cross by dragging the specific parent flies to the **Mating site** and by clicking on the **Mate** button.
- Mating results can be seen on the **Summary** tab and each resulting offspring can be viewed by clicking on the **Individual** tab.
- To start a new mating experiment, click on the **New experiment** button. Your current experiment will automatically be saved for you. To keep track of your experiments, you can rename the experiment within the **Saved experiments** window.

**Answer**

**3** From the results of some of your crosses, you suspect that some of the mutations are found in genes that are linked. You wish to determine which, if any, of these genes are linked. To help you with this task, the Mutant 5 strain has been provided. What is the genotype of Mutant 5? Verify this genotype by performing appropriate crosses.

*Note: Unlike the Mutant 1-4 strains, the Mutant 5 strain may contain more than one mutation and may not be true breeding.*

**Answer**

Blank area for student work.

**4** Identify which, if any, of these genes are linked and determine the map distances between linked genes. Describe the crosses that you did, their outcomes, and the conclusion that you drew from each of these crosses.

**Answer**

Blank area for student answer.

